


INDIAN FDA APPROVED COMPANY

AN ISO 13485 : 2016 COMPANY

CoreX

Posterior Cervical Occipito Reconstruction System



NEBULA®

A JOINTLESS JOINT

Since 1999

★ **Tit Corex Poly Screw**
(Fully / Partial Thread)

CAT NO. TC171



3.5
(Fully Thread)

4.0

CAT NO. TC172



4.0
(Partially Thread)

Screw **Size in mm**

3.5 10 to 34 (FULLY THREAD)

4.0 10 to 34 (FULLY THREAD)

4.0 20 to 34 (PARTIALY THREAD)

(In 2 mm variation)

★ **Tit Corex Occipital Screw**

CAT NO. TC173



4.0

4.5

Screw **Size in mm**

4.0 6 to 12 (2 mm Variation)

4.5



★ **Tit Corex Occipital Midline Plate**

CAT NO. TC174



Size


Small

Medium

Large

★ **Tit Corex Occipital Adjustable Plate**

CAT NO. TC175



Size

Small

Standard

★ **Tit Corex Connection Rod**

CAT NO. TC176



Rod

3.5

Size in mm

120, 240, 300

★ **Tit Corex Pre-curved Occipital Rod**

CAT NO. TC178



Rod

3.5

Size in mm

200

★ **Tit Corex Pre-curved Occipital Plate/Rod**
CAT NO. TC179



Rod Size in mm

3.5 200


★ **Tit Corex Dual Dia. Connection Rod**
CAT NO. TC180



Rod Size in mm

3.5 - 200
5.5 - 280 3.5 X 5.5 X 480

★ **Tit Corex Offset Domino**
CAT NO. TC183



Size

3.5 x 5.5


★ **Tit Corex Offset Triangle Domino**
CAT NO. TC184



Size

3.5 x 5.5

★ **Tit Corex Axial Domino**
CAT NO. TC185



Size

3.5 x 5.5

1. CAT NO. S3552
COREX DEPTH GAUGE



6. CAT NO. S3558
COREX ALIGNMENT TOOL Q.C.



2. CAT NO. S3553
COREX Q.C. DRILL BIT

For corex poly screw



For Occipital screw




7. CAT NO. S3560
COREX INNIE SCREW HOLDER



3. CAT NO. S3555
COREX Q.C. BONE TAP (CANCELLOUS TH.)

For corex poly screw



8. CAT NO. S3561
COREX ROD BENDER PAIR



4. CAT NO. S3556
COREX Q.C. BONE TAP (CORTICAL TH.)

For Occipital screw



9. CAT NO. S3562
COREX INNIE TIGHTNER Q.C.



5. CAT NO. S3557
COREX PEDICLE FEELER



10. CAT NO. S3563
COREX VERTEBRAL HOOK PROBE



11. CAT NO. S3564
COREX ROD TEMPLATE



13. CAT NO. S3571
COREX PEDICLE PROBE CU.



12. CAT NO. S3565
COREX ROD HOLDER FORCEP



14. CAT NO. S3572
COREX PEDICLE LOCALIZER



13. CAT NO. S3568
COREX POLY SCREW DRIVER



15. CAT NO. S3573
COREX ROD PUSHER



12. CAT NO. S3569
COREX ENTRY AWL



16. CAT NO. S3574
COREX DISTRACTOR CLAMP



13. CAT NO. S3570
COREX PEDICLE PROBE ST.




17. CAT NO. S3575
COREX ROD ROTATION FORCEP



18. CAT NO. S3577
COREX ROD STABILIZER (TORQUE)



19. CAT NO. S3578
COREX COMPRESSION CLAMP



20. CAT NO. S3579
COREX ROD REDUCER FORCEP




21. CAT NO. S3580
COREX QUICK COUPLING HANDLE ST.



22. CAT NO. S3581
COREX ROD CUTTER




23. CAT NO. S3582
COREX ROD BENDING PLIER






24. CAT NO. S3583
COREX QUICK COUPLE HANDLE T TYPE




25. CAT NO. 322.QC
SCREW DRIVER Q.C. 2.7MM FOR CROSS
LINK CONNECTOR



COREX INSTRUMENT BOX



COREX IMPLANT BOX





"Nebula" Narmada Park-3,
Opp. King's Height, VidhyaKunj Society Main Road,
Amin Marg, Rajkot-360005. Ph.: +91 281 2584832/33.



G-1921/3, Almighty Gate Main Road,
G.I.D.C Metoda-360021. Gujarat (INDIA)



Web : www.nebulasurgical.com



Email : info@nebulasurgical.com
export@nebulasurgical.com